



Med-Dyne Pediatric Disposable Pre-Wired Electrodes



Leading the Way in Pediatric Electrode Technology

2775 S. Floyd St.
Louisville, KY 40209

1-877-MED-DYNE **TOLL-FREE**
502-429-6759 **FAX**
sales@med-dyne.com **EMAIL**

www.med-dyne.com

Key Features

- **Patented Hydrogel Electrodes** are available in different shapes and sizes to meet your needs.
- **Color-Coded Lead Wires** made with Industry Standard Safety Sockets to ensure compatibility.
- **Pre-Wired Electrodes** make inventory levels easier to manage with one part for all of your ECG needs.

Convenient to Use

- 3 or 5 Lead sets in any length, shape, and size of choice.
- Electrodes are repositionable without need for added moisture.
- Each card contains the electrodes attached to the lead wires, all you need for the patient.

Color-Coded for Easy Connection

- Color-coded lead wires make cable connection easy.
- Strain relief device allows for more connection flexibility.

Comfortable to Wear


- Latex free.
- Breathable cloth for long term application.
- Flexible wires feel comfortable against patient's skin.
- More chest area is exposed to allow clinician to conduct other tests without disturbing the ECG monitoring.
- Our patented hydrogel layers provide firm adhesion while minimizing irritation to delicate newborn skin.

Ordering Information

Part #	Description	Quantity
G63-300	24" Prewired (BK,GR,WH) ECG Electrode Monitoring Neonate/ Pediatric Cloth Conductive Hydrogel, 7/8" x 1.25" Rectangle	300 Electrodes / 100 Cards
G63I-300	24" Prewired (BK,GR,WH) ECG Electrode Monitoring Neonate Cloth Conductive Hydrogel, 7/8" Square	300 Electrodes / 100 Cards
G63R-300	24" Prewired (BK,GR,WH) ECG Electrode Monitoring Neonate Cloth Conductive Hydrogel, 1.25" Diameter	300 Electrodes / 100 Cards
G63S40-300	40" Prewired (BK,GR,WH) ECG Electrode Monitoring Neonate Cloth Conductive Hydrogel, 1" Diameter	300 Electrodes / 100 Cards
G324S-300	24" Prewired (BK,GR,WH) ECG Electrode Monitoring Neonate Cloth Conductive Hydrogel, 1" Diameter	300 Electrodes / 100 Cards
G324R-300	24" Prewired (BK,GR,WH) ECG Electrode Monitoring Neonate Cloth Conductive Hydrogel, 1.25" Oval Diameter	300 Electrodes / 100 Cards
G324L-300	24" Prewired (BK,GR,WH) ECG Electrode Monitoring Neonate Cloth Conductive Hydrogel, 1.50" Diameter	300 Electrodes / 100 Cards



Part #	Description	Quantity
G524S-500	24" Prewired (BK,BR,GR,RD,WH) ECG Electrode Monitoring Neonate Cloth Conductive Hydrogel, 1" Diameter	500 Electrodes / 100 Cards
G524R-500	24" Prewired (BK,BR,GR,RD,WH) ECG Electrode Monitoring Neonate Cloth Conductive Hydrogel, 1.25" Oval Diameter	500 Electrodes / 100 Cards
G524L-500	24" Prewired (BK,BR,GR,RD,WH) ECG Electrode Monitoring Neonate Cloth Conductive Hydrogel, 1.50" Diameter	500 Electrodes / 100 Cards
G536S-500	36" Prewired (BK,BR,GR,RD,WH) ECG Electrode Monitoring Neonate Cloth Conductive Hydrogel, 1" Diameter	500 Electrodes / 100 Cards
G536L-500	36" Prewired (BK,BR,GR,RD,WH) ECG Electrode Monitoring Neonate Cloth Conductive Hydrogel, 1.50" Diameter	500 Electrodes / 100 Cards



Med-Dyne Value

Why use Pediatric Disposable Pre-Wired Electrodes?

- Significantly reduce your patient's risk of Healthcare-Associated Infections by using new electrodes and lead wires on each patient. 3 or 5 Lead sets in any length, shape, and size of choice.
- Americans get an estimated 1.7 million Healthcare-Associated Infections in the U.S. annually, according to the Centers for Disease Control and Prevention (CDC). Each year these infections:
 - ij Cause about 99,000 deaths.
 - ij Cost from \$35.7 to \$45 billion in additional healthcare costs annually.

Why choose Med-Dyne Pediatric Disposable Pre-Wired Electrodes?

- Med-Dyne's Pediatric Disposable Pre-Wired Electrodes are designed with a baby in mind.
 - ij Electrodes may be repositioned without using moisture.
 - ij Our patented hydrogel layers provide firm adhesion while minimizing irritation to delicate newborn skin.
 - ij Med-Dyne has the capability to custom design Pre-Wired Electrodes to fit your needs from color, length, to size.
- Med-Dyne always keeps your inventory in stock and ready to ship to reduce your on hand inventory costs.



2775 S. Floyd St.
Louisville, KY 40209

1-877-MED-DYNE **TOLL-FREE**
502-429-6759 **FAX**
sales@med-dyne.com **EMAIL**

www.med-dyne.com

Proudly made in the U.S.A.