

Clinical Monitoring Accessories



Gold Cup Electrode

Description:

Cloth gel electrode, 24" color-coded leads, black, white and green with a 1.5-mm (0.060") female connector.



Gold Cup Electrode

Lengths:

Standard lengths of 48", 60", 72", 96", and 120". Custom lengths and colors are available upon request.

Call for Pricing



Silver Cup Electrode

Lengths:

Standard lengths of 48", 60", 72", 96", and 120". Custom lengths and colors are available upon request.

Call for Pricing

Clinical Pediatric Pre-wired Electrodes

General Information:

Med-Dyne Clinical Pediatric pre-wired electrodes are designed with a baby in mind. Our patented hydrogel layers provide firm adhesion while minimizing irritation to delicate newborn skin. Clinical Pediatric electrodes may be repositioned without using moisture. All Clinical Pediatric electrodes have color-coded 24" lead wires with industry standard 0.060" (1.5mm) safety socket connector.

Features:

1. Small size
2. Tinned copper lead wire
3. Repositionable without using moisture
4. Breathable cloth for long term application
5. 0.060" (1.5mm) safety socket
6. Color-coded three or five leads available
7. Latex free to protect against latex sensitivity

G324S




Packaging:

3 per pouch

Dimensions: 0.9" Round

Call for Pricing

G524S



Description:

Includes Red, Black, White, Green and Brown leadwires

Packaging:

5 per pouch

Dimensions: 0.9" Round

Call for Pricing

G324R




Packaging:

3 per pouch

Dimensions: 1 1/4" Oval

Call for Pricing

G524R



Description:

Includes Red, Black, White, Green and Brown leadwires

Packaging:

5 per pouch

Dimensions: 1 1/4" Oval

Call for Pricing

G324L




Packaging:

3 per pouch

Dimensions: 1 1/2" Round

Call for Pricing

G524L



Description:

Includes Red, Black, White, Green and Brown leadwires

Packaging:

5 per pouch

Dimensions: 1 1/2" Round

Call for Pricing

Clinical Monitoring Accessories



DEL Series Snap Cloth Electrodes

Description:

Cloth snap, silver/silver chloride gel electrode

Uses:

Sleep, Infant Apnea, Pediatric / Adult applications

DEL3000



Packaging:

45 Electrodes Per Package.
10 Packages Per Case.

Dimensions: 7/8" x 1.25"

Call for Pricing

DEL3000A




Packaging:

45 Electrodes Per Package.
10 Packages Per Case.

Dimensions: 1.25" Square

Call for Pricing

DEL3000I



Packaging:

45 electrodes per package.
10 packages per case.


Dimensions: 7/8" square

Call for Pricing

Cables

Description:

We offer an extensive line of DIN style cables in 3.5 and 7 lead configurations with an over molded strain relief allowing for more flexibility thereby increasing durability. The yoke or patient end also is marked with color coding and letter designations.



Leadwires

Description:

DIN to snap lead wires have a 1.5 mm (0.060") over molded safety socket and at monitor end.



2775 South Floyd Street • Louisville, KY 40209-1817

1-877-MED-DYNE (633-3963) Phone

sales@med-dyne.com Email

502-429-6759 Fax